



Recombinant Human Early growth response protein 3 (EGR3)

| | |
|--------------------------|---|
| Product Code | CSB-BP007486HU |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q06889 |
| Form | Liquid or Lyophilized powder |
| Storage Buffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | ≥85% (SDS-PAGE) |
| Sequence | MTGKLAEKLPVTMSSLLNQLPDNLYPEEIPSA LNLFSGSSDSVVHYNQMATEN VMDIGLTNEKPNPELSYSGSFQPAPGNKTVTYLGKFAFDSPSNWCQDNIISLM SAGILGVPPASGALSTQTSTASMVQPPQGDVEAMYPALPPYSNCGDLYSEPV SFHDPQGNPGLAYSPQDYQSAKPA LDSNLFPMIPDYNLYHHPNDMGSIPEHK PFQGMPIRVNPPITPLETIKAFKDKQIHGPGFGLPQPPLTLKPIRPRKYPNRP SKTPLHERPHACPAEGCDRRFSRSDDELTRHLRIHTGHKPFQCRICMRSFSRSD HLTTHIRTHTGEKPFACEFCGRKFARSDERKRHAKIHLKQKEKKA EKG GAPS A SSAPPVSLAPVVTTCA |
| Source | Baculovirus |
| Target Names | EGR3 |
| Expression Region | 1-387 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.