



Recombinant Human Protocadherin-1 (PCDH1), partial

Product Code	CSB-YP017524HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q08174
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	TRVYKVPEEQPPNTLIGSLAADYGFPDVGHLYKLEVGAPYLRVDGKTGDIFTT ETSIDREGLRECQNQLPGDPCILEFEVSITDLVQNGSPRLLLEGQIEVQDINDNT PNFASPVITLAIPENTNIGSLFPIPLASDRDAGPNGVASVELQAGPEAQELFGLQ VAEDQEEKQPQLIVMGNLDLDRERWDSYDLTIKVDGGSPPRASSALLRVTVLDL NDNAPKFERPSYEAELSENSPIGHSVIQVKANDSDQGANAIEIYTFHQAPEVV RRLRLDRNTGLITVQGPVDREDLSTLRFVSLAKDRGTNPKSARAQVVTVKD MNDNAPTIEIRGIGLVTHQDGMANISEDVAEETAVALVQVSDRDEGENAAVTC VVAGDVPFQLRQASETGSDSKKKYFLQTTTPLDYEKVKDYTIEIVA VDSGNPPL SSTNSLKVQVVDVNDNAPVFTQSVTEVAFPENNKPGEVIAEITASDADSGSNA ELVYSLEPEPAAKGLFTISPETGEIQVKTSLDREQRESYELKVVAADRGSPSLQ GTATVLVNVLDCNDNDPKFMLSGYNFSVMENMPALSPVGMVTVIDGDKGENA QVQLSVEQDNGDFVIQNGTGTLSSLSFDREQQSTYTFQLKAVDGGVPPRSAY VGVTINVL DENDNAPYITAPSNTSHKLLTPQTRLGETVSQVAAEDFDSGVNAEL YSIAGGNPYGLFQIGSHSGAITLEKEIERRHHGLHRLVVKVSDRGGKPPRYGTA LVHLYVNETLANRTLLETLLGHSLDTPLDIDIAGDPEYERSKQRGN
Source	Yeast
Target Names	PCDH1
Expression Region	58-852
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.