



Recombinant Human Protein spire homolog 1 (SPIRE1)

Product Code	CSB-BP022588HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q08AE8
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AQAAGPAGGGEP RTEAVGGEGPREPGAAGGAAGGSRDALSLEEILRLYNQPI NEEQAWAVCYQCCGSLRAAARRRQPRHRVRSAAQIRVWRDGA VTLAPAADD AGEPPPVAGKLGYSQCMETEVIESLGI IYKALDYGLKENEERELSPPLEQLIDH MANTVEADGSNDEGYEAAEEGLGDEDEKRKISAIRSYRDVMKLCAAHLPTES DAPNHQAVCRALFAETMELHTFLTKIKSAKENLKKIQEMEKSDDESSTDLEELK NADWARFWVQVMRDLRNGVKLKKVQERQYNPLPIEYQLTPYEMLMDDIRCK RYTLRKVMVNGDIPRRLKKS AHEILD FIRSRPPLNPVSARKLKPTPPRPSLHE RILEEIK AERKLRPVSP EIRRSRLAMRPLSMSYSFDLSDVTTPESTKNLVESS MVNGGLTSQTKENGLSTSQQVPAQRKLLRAPTLAELDSS ESEEEETLHKSTSS SSVSPSFPEEPVLEAVSTRKKPKFLPISSTPQPERRQPPQRRHSIEKETPTNV RQFLPPSRQSSRSLEEF CYPVECLALTVEEVMHIRQVLVKA ELEKYQQYKDIYT ALKKGKLCFCCRTRRF SFFTWSYTCQFCKRPVCSQCCKMRLPSKPYSTLP IF SLGPSALQRGESSMRSEKPSTAHHRPLRSIARFSSKSKSMDKSDEELQFPKEL MEDWSTMEVCVDCKKFISEI ISSSRSLVLANKRARLKRKTQSFYMS SPGPSE YCPSERTISEI
Source	Baculovirus
Target Names	SPIRE1
Expression Region	2-756
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.