



Recombinant Human E3 ubiquitin-protein ligase DTX4 (DTX4)

Product Code	CSB-MP764784HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y2E6
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLLASAVVVWEWLNEHGRWRPYSPAVSHHIEAVVRAGPRAGGSVVLGQVDS RLAPYIIDLQSMNQFRQDTGTLRPVRRNYYDPSSAPGKGVVWEWENDNGSW TPYDMEVGITIQHAYEKQHPWIDLTSIGFSYVIDFNTMGQINRQTQRQRRVRRR LDLIYPMVTGTLPKAQSWPVSPGPATSPPMSPCSCPQCVLVMSVKAADVNGS TGPLQLPVTRKNMPPPGVVKLPLPGSGAKPLDSTGTIRGPLKTAPSQVIRRQ ASSMPTGTTMGSPASPPGPNSTGRVALATLNRTNLQRLAIAQSRVLIASGVP TVPVKNLNGSSPVNPALAGITGILMSAAGLPVCLTRPPKLVLHPPVSKSEIKSI PGVSNTSRKTTKKQAKKGKTPEEVLLKYLQKVRHPPDEDCTICMERLTAPSGY KGPQPTVKPDLVGKLSRCGHVYHIYCLVAMYNNGNKDGLQCPTCKTIYGVK TGTQPPGKMEYHLIPHSPLGHPDCKTIRIISIPPGIQGPEHPNPGKSFSARGFP RHCYLPDSEKGRKVLKLLVAVDRRLIFAIGTSSTTGESDVIWNEVHHKTEFG SNLTGHGYPDANYLDNVLAELAAQGISDSTAQEKD
Source	Mammalian cell
Target Names	DTX4
Expression Region	1-619
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.