



# Recombinant Human Tubulin monoglycylase TTLL3 (TTLL3)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP2332HU-B   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | Q9Y4R7   |
| <b>Form</b>              | Liquid or Lyophilized powder   |
| <b>Storage Buffer</b>    | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Homo sapiens (Human)   |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MNRLRNAKIYVERAVKQKKIFTIQGCYPVIRCLLRRRGWVEKMMVHRSGPTLL<br>PPQKDL DSSAMGDS DTT EDEDEDEDEEFQPSQLFDFD DLLKFDDL DGTALM<br>VGLCLNLRNLPWFDEVDANSFFPRCYCLGAEDDKAFIEDFWLTAARNVLKLV<br>VKSEWKSYP IQAVEEEASGDKQPKKQEKNPVLSPEFVDEALCACEEYLSNLA<br>HMDIDKDLEAPLYLTPEGWSLFLQRY YQVVHEGAELRHLDTQVQRCE DILQQL<br>QAVVPQIDMEGDRNIWIVKPGAKSRGRGIMCMDHLEMLKLVNGNPVVMKDG<br>KWVVQKYIERPLLIFGTKFDLRQWFLVTDWNPLTVWFYRDSYIRFSTQPFSLK<br>NLDNSVHLCNNSIQKHL ENSCHRHPLLPPDNMWSSQRFQAH LQEMGAPNAW<br>STIIVPGMKDAVIHALQTSQDTVQCRKASFELYGADVFVGFEDFQPW LIEINASP<br>TMAPSTAVTARLCAGVQADTLRVVIDRMLDRNC DTGAFELIYKQPAVEVPQYV<br>GIRLLVEGFTIKKPMAMCHRRMGVRPAVPLLTQRGSGEARHHFPSLHTKAQLP<br>SPHVL RHQQVLR RQH SKLVG TKALSTTGKALRTLPTAKVFISLPPNLDFKVP<br>SILKPRKAPALLCLRGPQLEVPCCLCPLKSEQFLAPVGRSRPKANSRPDCDKP<br>RAEACPMKRLSPLKPLPLVGTFQRRRGLGDMKLGKPLLRFP TALVLDPTPNKK<br>KQVKYLGLDSIAVGGSRVDGAR PCTPGSTAR A |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | TTLL3  |
| <b>Expression Region</b> | 1-772  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a   |



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.