



Recombinant Human Nuclear protein AMMECR1 (AMMECR1)

Product Code	CSB-EP2255HU-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y4X0
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAAGCCGVKKQKLSSSPSGSGGGGGASSSSHCSGESQCRAGELGLGGAG TRLNGLGGLTGGGSGSGCTLSPPQGC GG GGGIALSPPPSCGVGTLLSTPAA ATSSSPSSSSAASSSSPGSRKMVSAEMCCFCFDVLYCHLYGYQQPRTPRFT NEPYPLFVTWKIGRDKRLRGCI GTF SAMNLHSGLREYTLTSALKDSRFP PMTR DELPRLFCSV SLLTNFEDVCDYLDWEVGVHGIRIEFINEKGSKRTATYLPEVAK EQGWDHIQTIDSLLRKGGYKAPITNEFRKTIKLTRYRSEKMTLSYAEYLAHRQH HHFQNGIGHPLPPYNHYS
Source	E.coli
Target Names	AMMECR1
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.