



Recombinant *Thermotoga maritima* 5-methylthioadenosine/S-adenosylhomocysteine deaminase (mtaD)


Product Code	CSB-EP895804TNJ
Abbreviation	Recombinant <i>Thermotoga maritima</i> mtaD protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9X034
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant <i>Thermotoga maritima</i> 5-methylthioadenosine/S-adenosylhomocysteine deaminase (mtaD)
Immunogen Species	<i>Thermotoga maritima</i> (strain ATCC 43589 / DSM 3109 / JCM 10099 / NBRC 100826 / MSB8)
Sensitivity	Not Test
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MIIGNCLILKDFSSEPFWGAVEIENGTIKRVLQGEVKVDLDLSGKLVMPALFNTH THAPMTLLRGVAEDLSFEEWLFVSKVLPIDRLTEKMAYYGTLAQMEMARHGIA GFVDMYFHEEWIAKAVRDFGMRALLTRGLVDSNGDDGGRLEENLKLYNEWN GFEGRIFVGFPHSPYLCSEEYLRVFDATAKSLNAPVTIHLHYETSKEEYDLEDIL NIGLKEVKTIAAHCVHLPERYFGVLKDIPFFVSHNPASNKLKLGNGIAPVQRMIEH GMKVTLGTDGAASNNSLNLFFEMRLASLLQKAQNPRNLDVNTCLKMVTYDGA QAMGFKSGKIEEGWNADLVVIDLDLPEMFVPVQNIKNHLVHAFSGEVFATMVAG KWIYFDGEYPTIDSEEVKRELARIEKELYSS
Research Area	Others
Source	E.coli
Target Names	mtaD
Expression Region	1-406aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	52.7 kDa



Protein Length

Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.