





# Recombinant *Macaca fascicularis* Interleukin-2 (IL2) (Active)


<b>Product Code</b>	CSB-MP640916MOVd9
<b>Abbreviation</b>	Recombinant Cynomolgus monkey IL2 protein (Active)
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q29615
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 μm sterile filtered PBS, 6% Trehalose, pH 7.4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Macaca fascicularis</i> (Crab-eating macaque) ( <i>Cynomolgus monkey</i> )
<b>Biological Activity</b>	① Measured by its binding ability in a functional ELISA. Immobilized <i>Macaca fascicularis</i> IL2RA (CSB-MP5104MOV) at 2 μg/mL can bind <i>Macaca fascicularis</i> IL2. The EC50 is 88.60-116.8 ng/mL. ② <i>Macaca fascicularis</i> IL2 captured on Protein A Chip can bind Recombinant <i>Macaca fascicularis</i> IL2RA (CSB-MP5104MOV) with an affinity constant of 31.3 nM as detected by MetaSPR Assay (WeSPRTM 200).
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	APTSSSTKKTQLQLEHLLLDLQMILNGINNYKNPKLTRMLTFKFYMPKKATELR HLQCLEEELKPLEEVLNLAQSKSFHLRDTKDLISNINVIVLELKGSETTLMCEYA DETATIVEFLNRWITFCQSIISTLT
<b>Source</b>	Mammalian cell
<b>Target Names</b>	IL2
<b>Expression Region</b>	21-154aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal hFc1-tagged
<b>Mol. Weight</b>	44.4 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



**Activity**  
Measured by its binding ability in a functional ELISA. Immobilized *Macaca fascicularis* IL2RA (CSB-MP5104MOV) at 2 µg/ml can bind *Macaca fascicularis* IL2. The EC<sub>50</sub> is 88.60-116.8 ng/mL.



**Activity**  
*Macaca fascicularis* IL2 captured on Protein A Chip can bind Recombinant *Macaca fascicularis* IL2RA (CSB-MP5104MOV) with an affinity constant of 31.3 nM as detected by MetaSPR Assay (WeSPRTM 200).

### Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.