




Recombinant Pongo abelii Glial fibrillary acidic protein (GFAP)

Product Code	CSB-EP009369PYXe1
Abbreviation	Recombinant Pongo abelii GFAP protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5RA72
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Pongo abelii Glial fibrillary acidic protein (GFAP)
Immunogen Species	Pongo abelii (Sumatran orangutan) (Pongo pygmaeus abelii)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	MERRRITSAARRSYISSGEMMVGGLAPGRRLGPGTRLRLSLARMPPPLPTRVDF SLAGALNAGFKETRASERAEMMELNDRFASYIEKVRFLQQNKALAAELNQLR AKEPTKLADVYQAELELRLRLDQLTANSARLEVERDNLAQDLGTVRQKLQDE TNLRLEAENNLAAAYRQEADATLARLDLERKTESLEEEIRFLRKIHEEEVRELQ EQLARQQVHVELDMAKPDLTAAALKEIRTQYEAMASSNMHEAEEWYRSKFADL TDAAARNAELLRQAKHEANDYRRQLQSLTCDLESRLRGNTNESLERQMREQEER HAREAASYQEALARLEEEEGQSLKDEMARHLQEYQDLLNVKLALDIEIATYRKL EGEENRITIPVQTFSNLQIRETSLDTKSVSEGHKRNIVVKTVMRDGEVIKESK QEHKDVM
Research Area	Neuroscience
Source	E.coli
Target Names	GFAP
Expression Region	1-432aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	50.0 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.