




Recombinant Cat Prolactin (PRL)

Product Code	CSB-EP018724CA
Abbreviation	Recombinant Cat PRL protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46403
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Cat Prolactin (PRL)
Immunogen Species	Felis catus (Cat)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	LPICPSGAVNCQVSLRDLFDRAVILSHYIHNLSSSEMNFNEFDKRYAQGRGFITKAINSCHTSSLPTPEDKEQAQQIHHEDLLNVILRVLRSWNDPLYHLVTEVRGLHEAPDAILSRAIEIEEQNRRLLEGMEKIVHQVHPGVRENEVYSVWSGLPSLQMADEDSRLFAFYNLLHCLRRDSHKIDSYLKLLKCRIVYDSNC
Research Area	Signal Transduction
Source	E.coli
Target Names	PRL
Expression Region	31-229aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	29.9 kDa
Protein Length	Full Length of Mature Protein
Image	 <p>(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.</p>
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$. The shelf life of lyophilized form is 12 months at $-20^{\circ}\text{C}/-80^{\circ}\text{C}$.