



Recombinant Human Protein-glutamine gamma-glutamyltransferase 2 (TGM2) (Y443G)


Product Code	CSB-EP023462HU(A4)(M2)
Abbreviation	Recombinant Human TGM2 protein (Y443G)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21980
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human Protein-glutamine gamma-glutamyltransferase 2 (TGM2) (Y443G)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MAEELVLERCDLELETNGRDHHTADLCREKLVVRRGQPFWLTLHFEGRNYEA SVDSLTFVVTGPAPSQEAGTKARFPLRDAVEEGDWTATVVDQQDCTLSLQL TTPANAPIGLYRLSLEASTGYQGSSFVLGHFILLFNAWCPADAVYLDSEEERQE YVLTQQGFIYQGSAKFIKNIPWNFGQFEDGILDICLILLDVNPKFLKNAGRDCSR RSSPVYVGRVVS GMVNCNDDQGVLLGRWDNNYGDGVSPMSWIGSVDILRR WKNHGCQRVKYGQCWVFAAVACTVLRCLGIPTRVVTNYNSAHDQNSNLLIEY FRNEFGEIQGDKSEMIWNFHCWVESWMTRPDLQPGYEGWQALDPTPQEKSE GTYCCGPVPVRAIKEGDLSTKYDAPFVFAEVNADVVDWIQQDDGSVHKSINRS LIVGLKISTKSVGRDEREDITHGKYPEGSSEEREAFTANHLNKLAEKEETGM AMRIRVGQSMNMGSDFDVFAHITNNTAEYVCRLLLCARTVSYNGILGPECGT KYLNLNLEPFSEKSVPLCILEYKYRDCLTESNLIKVRALLVEPVINSYLLAERDL YLENPEIKIRILGEPKQKRKLVAEVSLQNPLVALEGCTFTVEGAGLTEEQKTVE IPDPVEAGEEVKVRMDLLPLHMGLHKLVVNFESDKLKAVKGFNRNVIIGPA
Research Area	Immunology
Source	E.coli
Target Names	TGM2
Expression Region	1-687aa(Y443G)
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged



Mol. Weight 84.1 kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.