




Recombinant Mouse Cadherin-16 (Cdh16), partial

Product Code	CSB-MP005041MO
Abbreviation	Recombinant Mouse Cdh16 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O88338
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Mouse Cadherin-16 (Cdh16), partial
Immunogen Species	Mus musculus (Mouse)
Sensitivity	Not Test
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	QPAESLHTEVPENYGGNFPFYILKLPLPLGRDEGHIVLSGDSNTADQNTFAVDT DSGFLVATRTLDRREEKAEYQLQVTLESEDGRILWGPQLVTVHVKDENDQVPQ FSQAIYRAQLSQGTRPGVPFLFLEASDGDAPGTANSDLRFHILSQSPPQPLPD MFQLDPHLGALALSPSGSTSLDHALEETYQLLVQVKMDGDQPSGHQAIATVEI SIVENSWAPLEPVHLAENLKVVYPHSIAQVHWSSGGDVHYQLESQPPGPFVDV TEGMLHVTMELDREAQAEYQLQVRAQNSHGEDIAEPLQLQVVMDENDNAP VCSPHDPTVNIPELSPPGTEIARLSAEDLDAPGSPNSHIVYQLLSPEPEEGAEN KAFELDPTSGSVTLGTAPLHAGQSILLQVLAVDLAGSEGLSSTCEVTVMVTDV NNHAPEFINSQIGPVTLPEDVKPGALVATLMATDADLEPAFRLMDFAIEEGDPE GIFDLSEWPDSDHVQLRLRKNLSYEAAPDHKVVVVVSNIEELVGP GP GPAATA TVTILVERVVAPLKLDDQESYETSIPVSTPAGSLLLTIQPSDPMSTLRFSLVNDS EGWLCIKEVSGEVHTAQSLQGAQPGDTYTVLVEAQDTPGLSTSATVVIHFL KASPVPALTL SAGPSRHLCTPRQDYGVVSGVSEDPDLANRNGPYSFALGPN PTVQRDWRLQPLNDSHAYLTALHWWEPGEYMPV VVHH DTHMWQLQVKVI VCRCNVEGQCMRKVGRMKGMPTKLSA
Research Area	Signal Transduction
Source	Mammalian cell
Target Names	Cdh16
Expression Region	22-788aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	84.9 kDa


Protein Length

Partial

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.