




Recombinant Human V-type proton ATPase subunit B,brain isoform (ATP6V1B2)

Product Code	CSB-EP002398HUc7
Abbreviation	Recombinant Human ATP6V1B2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21281
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Human V-type proton ATPase subunit B,brain isoform (ATP6V1B2)
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MALRAMRGIVNGAAPELPVPTGGPAVGAREQALAVSRNYLSQPRLTYKTVSG VNGPLVILDHVKFPARYAEIVHLTLPDGTKRSGQVLEVSGSKAVVQVFEGTSGID AKKTSCEFTGDILRTPVSEDM LGRVFN GSGKPIDRGPV VLAEDFLDIMGQPINP QCRIYPEEMIQTGISAIDGMNSIARGQKIPIFSAAGLPHNEIAAQICRQAGLVKKS KD VVDYSEENFAIVFAAMGVN METARFFKSDFEENGSM DN VCLFLN LANDPTI ERIITPRLALTTAEFLAYQCEKHV LVILTDMSSYAEALREVSAA REEVPGRRGFP GYMYTDLATIYERAGRVEGRNGSITQIPILTMPND DITHPIPDLTGYITEGQIYVD RQLHNRQIYPPINVLPSLSRLMKSAIGEGMTRKDHADVSNQLYACYAIGKD VQ AMKAVVGEEALTSDDL LYLEFLQKFERNFIAQGPYENRTVFETLDIGWQLLRIF PKEMLK RIPQSTLSEFYPRDS AKH
Research Area	Metabolism
Source	E.coli
Target Names	ATP6V1B2
Expression Region	1-511aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-tagged
Mol. Weight	63.4 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.