



Recombinant Human Dipeptidyl peptidase 4 (DPP4), partial

Product Code	CSB-MP007139HU1d7
Abbreviation	Recombinant Human DPP4 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27487
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE. ≥ 95% as determined by SEC-HPLC.
Sequence	NKGTDDATADSRKTYTLTDYLNKNTYRLKLYSLRWISDHEYLYKQENNILVFNAE YGNSSVFLENSTFDEFGHSINDYSISPDGQFILLEYNVVKQWRHSYASYDIYD LNKRQLITEERIPNNTQWVTWSPVGHKLAYVWNNDIYVKIEPNLPSYRITWTGK EDIIYNGITDWVYEEVFSAYSALWWSPNGTFLAYAQFNDTEVPLIEYSFYSD SLQYPKTVRVPYPKAGAVNPTVKFFVNTDSSLSSVTNATSIQITAPASMLIGDH YLCDVTWATQERISLQWLRRIQNYVMIDICDYDESSGRWNCLVARQHIEMST TGWVGRFRPSEPHFTLDGNSFYKIISNEEGYRHICYFQIDKKDCTFITKGTWEVI GIEALTSYLYYISNEYKGMPPGGRNLYKIQLSDYTKVTCLSCELNPERCQYYSV SFSKEAKYYQLRCSGPGPLPLYTLHSSVNDKGLRVLEDNSALDKMLQNVQMPS KKLDFIILNETKFWYQMILPPHFDKSKKYPLLLDVYAGPCSQKADTVFRLNWAT YLASTENIIVASFDGRGSGYQGDKIMHAINRRLGTFEVEDQIEAARQFSKMGFV DNKRIAIWGSYGGYVTSMLVLSGSGVFKCGI AVAPVSRWEY YDSVYTERYM GLPTPEDNLDHYRNSTVMSRAENFKQVEYLLIHGTADDNVHFQQSAQISKALV DVGVD FQAMWYTDDEDHGIASSTAHQH IYTHMSHF I KQC FSLP
Research Area	Immunology
Source	Mammalian cell
Target Names	DPP4
Expression Region	29-766aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.




Tag Info C-terminal 10xHis-tagged


Mol. Weight 86.7 kDa

Protein Length Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of DPP4 was greater than 95% as determined by SEC-HPLC

Reconstitution We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.