




Recombinant Human Neogenin (NEO1), partial

| | |
|--------------------------|---|
| Product Code | CSB-MP852917HU |
| Abbreviation | Recombinant Human NEO1 protein, partial |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q92859 |
| Form | Liquid or Lyophilized powder |
| Storage Buffer | If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0. |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Sensitivity | Not Test |
| Purity | ≥ 95% as determined by SDS-PAGE. |
| Sequence | <p>ARSGSAPQSPGASIRTFYFLVEPVDTLVSRGSSVILNCSAYSEPSPKIEWK KDGTFLNLVSDRRQLLPDGSFLISNVVHSHKHNKPDEGYQCVAIVESLGTIIS RTAKLIVAGLPRFTSQPEPSSVYAGNNAILNCEVNADLVPFVRWEQNRQPLLL DDRVIKLPSPGMLVISNATEGDGGLYRCVVESSGGPPKYSDEVELKVLDPPEVIS DLVFLKQPSPLVRVIGQDVVLPVASGLPTPTIKWMKNEEALDTESSERLVLLA GGSLEISDVTEDDAGTYFCIADNGNETIEAQAELTVQAQPEFLKQPTNIYAHES MDIVFECEVTGKPTPTVKWVKNMGDMVIPSDFKIVKEHNLQVLGLVKSDEGFY QCIAENDVGNAQAGAQLIILEHAPATTGPLPSAPRDVVASLVSTRFIKLTWRTP ASDPHGDNLTYSVFYTKEGIARERVENTSHPGEMQVTIQNLMPATVYIFRVMA QNKHSGGESSAPLRVETQPEVQLPGPAPNLRAYAASPTSITVTWETPVSGNG EIQNYKLYYMEKGTDKQDQDVSSHSYTINGLKKYTEYSFRVVAYNKHGPGVS TPDVAVRTLSDVPSAAPQNLVSRNSKSMIHWQPPAPATQNGQITGYKIRYR KASRKSDVTETLVSGTQLSQLIEGLDRGTEYNFRVAALTINGTGPATDWLSAET FESDLDETRVPEVPSSLHVRPLVTSIVVSWTPPENQIVVRGYAIGYGIGSPHA QTIKVDYKQRYTYIENLDPSSHVITLKAFNNVGEIPLYESAIVTRPHTDTSEVD LFINAPYTPVPDPTPMMPPVGVQASILSHDTIRITWADNSLPKHQKITDSRYYT VRWKTNIPANTKYKNANATTLVSYLVTLKPNLTYEFSVMVTKGRRSSTWSMTA HGTTFELVPTSPPKDVTVVSKEGKPKTIIVNWQPPSEANGKITGYIIYYSTDVNA EIHDWVIEPVVGNRLTHQIQELTLTPYFYFKIQARNKGMGPMSEAVQFRTPK ADSSDKMPNDQASGSGGKGSRLPDLGSDYKPPMSGNSPHGSPTSPLDSN ML</p> |
| Research Area | Neuroscience |
| Source | Mammalian cell |
| Target Names | NEO1 |



| | |
|--------------------------|---|
| Expression Region | 34-1105aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | C-terminal hFc1-tagged |
| Mol. Weight | 146.5 kDa |
| Protein Length | Partial |

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.