



Phospho-IMPDH2 (S416) Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA011701XA01phHU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P12268 |
| Immunogen | Peptide sequence around site of S (416) derived from Human IMDH2 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, WB |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | IMPDH2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |