




PDS5B Antibody

Product Code	CSB-PA7047XA21DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	F4I735
Immunogen	Recombinant Arabidopsis thaliana (Mouse-ear cress)PDS5B protein (1107-1122aa)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana (Mouse-ear cress)
Tested Applications	ELISA, DB; Recommended dilution: DB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Target Names	PDS5B

Image



Dot blot analysis of peptides corresponding to the immunogen was probed with PDS5B antibody at a dilution of 1:1000. Blocking/Dilution buffer: 4% nonfat dry milk(NFDM)/PBS. Secondary antibody?Goat polyclonal to rabbit IgG at 1/40000 dilution Exposure time: 14s.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.