




NEP Antibody

Product Code	CSB-PA7391XA01IPN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D3J181
Immunogen	Recombinant Influenza A virus (A/reassortant/NIBRG-14(Viet Nam/1194/2004 x Puerto Rico/8/1934)(H5N1))NEP protein (6-116aa)
Raised In	Rabbit
Species Reactivity	Influenza A virus (A/reassortant/NIBRG-14(Viet Nam/1194/2004 x Puerto Rico/8/1934)(H5N1))
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Influenza A virus (A/reassortant/NIBRG-14(Viet Nam/1194/2004 x Puerto Rico/8/1934)(H5N1))
Target Names	NEP

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: MUS81 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 36.5 kDa

Observed band size: 36.5 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.