




Venom allergen 5.01 Antibody

Product Code	CSB-PA330485XA01VBP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35781
Immunogen	Recombinant Vespa crabro (European hornet)Venom allergen 5.01 protein (1-202aa)
Raised In	Rabbit
Species Reactivity	Vespa crabro (European hornet)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vespa crabro (European hornet)
Target Names	Venom allergen 5.01

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes:Venom allergen 5.01 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/20000 dilution

Predicted band size: 29.6 kDa

Observed band size: 29.6 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.