




TNFRSF25 Recombinant Monoclonal Antibody

Product Code	CSB-RA839000MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93038
Immunogen	Recombinant Human TNFRSF25 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG4
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFRSF25
Clone No.	hz0072

Image



The Binding Activity of Human TNFRSF25 with Anti-TNFRSF25 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TNFRSF25(CSB-MP839000HUn6) at 2 µg/mL can bind Anti-TNFRSF25 recombinant antibody. The EC₅₀ is 4.485-5.277 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.