




CD8A Recombinant Monoclonal Antibody

Product Code	CSB-RA004966MA3HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01732
Immunogen	Recombinant Human CD8A protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG2
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CD8A
Clone No.	9F1

Image



The Binding Activity of Human CD8A with Anti-CD8A Recombinant antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CD8A (CSB-MP004966HUd7) at 5 µg/mL can bind Anti-CD8A recombinant antibody. The EC₅₀ is 11.49-20.78 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.