



VEGFA (Brolucizumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA025833MB8HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15692
Immunogen	Recombinant Human VEGFA protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Human IgG1 kappa
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	VEGFA
Accession NO.	6C3

Image



The Binding Activity of Human VEGFA with Anti-VEGFA recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human VEGFA (CSB-YP025833HU) at 2 µg/mL can bind Anti-VEGFA recombinant antibody. The EC₅₀ is 26.68 - 30.32 ng/mL.