




TNFRSF17 Recombinant Monoclonal Antibody

Product Code	CSB-RA023974MA4HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02223
Immunogen	Recombinant Human TNFRSF17 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TNFRSF17
Clone No.	4D7E5

Image



The Binding Activity of Human TNFRSF17 with Anti-TNFRSF17 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TNFRSF17 protein (CSB-MP023974HU1d7) at 2 µg/mL can bind Anti-TNFRSF17 recombinant antibody. The EC₅₀ is 14.97-19.19 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.