




CSF1R Recombinant Monoclonal Antibody

Product Code	CSB-RA006044MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07333
Immunogen	Recombinant Human CSF1R protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG4
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CSF1R
Clone No.	18E11

Image



The Binding Activity of Human CSF1R with Anti-CSF1R recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CSF1R (CSB-MP006044HUd7) at 2 µg/mL can bind Anti-CSF1R recombinant antibody. The EC₅₀ is 1.800-2.218 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.