




FCGR2A Recombinant Monoclonal Antibody

Product Code	CSB-RA008540MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12318
Immunogen	Recombinant Human FCGR2A protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FCGR2A
Clone No.	5H2

Image



The Binding Activity of Human FCGR2A with Anti-FCGR2A Recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human FCGR2A (CSB-MP008540HU1) at 2 µg/mL can bind Anti-FCGR2A recombinant antibody. The EC₅₀ is 24.44-32.57 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.