




LGR5 Recombinant Monoclonal Antibody

Product Code	CSB-RA012906MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75473
Immunogen	Recombinant Human LGR5 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG2
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LGR5
Clone No.	7F1

Image



The Binding Activity of Human LGR5 with Anti-LGR5 Recombinant antibody
 Activity: Measured by its binding ability in a functional ELISA. Immobilized Human LGR5 (CSB-MP012906HU(A4)) at 10 µg/mL can bind Anti-LGR5 recombinant antibody. The EC₅₀ is 1.690-1.867 ng/mL. The VLPs (CSB-MP3838) is negative control.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.