




IL1RL1 Recombinant Monoclonal Antibody

Product Code	CSB-RA011626MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01638
Immunogen	Recombinant Human IL1RL1 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG2
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL1RL1
Clone No.	8G2

Image



The Binding Activity of Human IL1RL1 with Anti-IL1RL1 recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human IL1RL1(CSB-MP011626HU) at 2 µg/mL can bind Anti-IL1RL1 recombinant antibody. The EC₅₀ is 66.366-7.698 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.