





CRLF2 Recombinant Monoclonal Antibody

Product Code	CSB-RA005983MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HC73
Immunogen	Recombinant Human CRLF2 protein
Species Reactivity	Human, Macaca fascicularis
Tested Applications	ELISA, WB; Recommended dilution: ELISA:1:5000-1:50000, WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CRLF2
Clone No.	10A4


Image



The Binding Activity of Macaca fascicularis CRLF2 with Anti-CRLF2 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Macaca fascicularis CRLF2(CSB-MP5357MOVd7) at 2 µg/mL can bind Anti-CRLF2 recombinant antibody. The EC₅₀ is 2.079-2.450 ng/mL.



The Binding Activity of Human CRLF2 with Anti-CRLF2 recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CRLF2(CSB-MP005983HU-B) at 2 µg/mL on streptavidin coated plates can bind Anti-CRLF2 recombinant antibody. The EC₅₀ is 3.631-7.668 ng/mL.



Western Blot

Positive WB detected in: Jurkat whole cell lysate(30μg), THP-1 whole cell lysate(30μg), A-431 whole cell lysate(30μg), SH-SY5Y whole cell lysate(30μg)

All lanes: CRLF2 antibody at 1:1000


Secondary

Goat polyclonal to human IgG at 1/40000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa

Exposure time: 2min



The Binding Activity of Human CRLF2 with Anti-CRLF2 recombinant antibody

Activity: Measured by its binding ability in a functional ELISA. Immobilized Human CRLF2(CSB-MP005983HU). at 2 μg/mL can bind Anti-CRLF2 recombinant antibody The EC₅₀ is 1.630-1.842 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.