




SLC2A1 Recombinant Monoclonal Antibody

Product Code	CSB-RA021546MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11166
Immunogen	Recombinant Human SLC2A1 protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SLC2A1
Clone No.	18A8

Image



The Binding Activity of Human SLC2A1 with Anti-SLC2A1 Recombinant antibody Activity: Measured by its binding ability in a functional ELISA. Immobilized Human SLC2A1(CSB-MP021546HU(A4)d7) at 10 µg/mL can bind Anti-SLC2A1 recombinant antibody. The EC₅₀ is 5.023-6.275 ng/mL. The VLPs (CSB-MP3838) is negative control.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.