




TSHR Recombinant Monoclonal Antibody

Product Code	CSB-RA025131MA2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16473
Immunogen	Recombinant Human TSHR protein
Species Reactivity	Human
Tested Applications	ELISA; Recommended dilution: ELISA:1:5000-1:50000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	TSHR
Clone No.	K1-70

Image



The Binding Activity of Human TSHR with Anti-TSHR Recombinant antibody
Activity: Measured by its binding ability in a functional ELISA. Immobilized Human TSHR (CSB-MP025131HU) at 10 µg/mL can bind Anti-TSHR recombinant antibody. The EC₅₀ is 37.69-53.08 ng/mL.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.