



IL6R (Sapelizumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA011665MB2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08887
Immunogen	Recombinant Human IL-6R protein
Species Reactivity	Human
Tested Applications	ELISA, FACS
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	IL-6R
Accession NO.	5E9