




FGF13 Recombinant Monoclonal Antibody

Product Code	CSB-RA008619MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92913
Immunogen	Recombinant Human FGF13 protein
Species Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FGF13
Clone No.	14B1

Image



Western Blot

Positive WB detected in: Mouse Brain tissue lysate(30µg), Rat brain tissue lysate(30µg)
All lanes: FGF13 antibody at 1:1000

Secondary

Goat polyclonal to human IgG at 1/40000 dilution
Predicted band size: 28 kDa
Observed band size: 28 kDa
Exposure time: 2min

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.