




LRRTM4 Recombinant Monoclonal Antibody

Product Code	CSB-RA774813MA1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86VH4
Immunogen	Recombinant Human LRRTM4 protein
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Affinity-chromatography
Isotype	hIgG1
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	LRRTM4
Clone No.	21D2

Image




Western Blot

Positive WB detected in: HEK293 whole cell lysate(30µg),MCF7 whole cell lysate(30µg)
All lanes: LRRTM4 antibody at 1:1000

Secondary

Goat polyclonal to human IgG at 1/40000 dilution
Predicted band size: 67 kDa
Observed band size: 80 kDa
Exposure time: 120s



IHC image of CSB-RA774813MA1HU diluted at 1:50 and staining in paraffin-embedded human brain tissue performed on a Leica BondTM system. After dewaxing and hydration, antigen retrieval was mediated by high pressure in a citrate buffer (pH 6.0). Section was blocked with 10% normal goat serum 30min at RT. Then primary antibody (1% BSA) was incubated at 4°C overnight. The primary is detected by a Anti-Human IgG, Fcy Fragment Specific labeled by HRP and visualized using 0.05% DAB.



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.