




Recombinant Aspergillus niger Probable endo-1,4-beta-xylanase B (xlnB)

Product Code	CSB-YP383243AVEf3
Abbreviation	Recombinant Aspergillus niger xlnB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A2Q7I0
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Aspergillus niger (strain ATCC MYA-4892 / CBS 513.88 / FGSC A1513)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	VPHDSVAQRSDALHMLSERSTPSSTGENNGFYYSFWTDGGGDVTYTNGDAG AYTVEWSNVGNFVGGKGWNPQSAQDITYSGTFTPSGNGYLSVYGWTTDPLIE YYIVESYGDYNPGSGGTYKGTVTSDGSVYDIYTATRTNAASIQTATFTQYWS VRQNKRVGGTVTTSNHFNAWAKLGMNLGTHNYQIVATEGYQSSGSSSITVQ
Research Area	Others
Source	Aspergillus niger (strain ATCC MYA-4892 / CBS 513.88 / FGSC A1513)
Target Names	xlnB
Expression Region	19-225aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 6xHis-Myc-tagged
Mol. Weight	25.9 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.