



Recombinant Human Maltase-glucoamylase (MGAM), partial

Product Code	CSB-MP013772HU(A4)
Abbreviation	Recombinant Human MGAM protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O43451
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	SAECPVVNELERINCIPDQPPTKATCDQRGCCWNPQGAVSVPWCYYSKNHS YHVEGNLVNTNAGFTARLKNLPSSPVFGSNVDNVLTAEYQTSNRFHFKLTDQ TNNRFEVPHEHVQSFSGNAAASLTYQVEISRQPFSSIKVTRRSNNRVLFDSSIGP LLFADQFLQLSTRLPSTNVYGLGEHVHQYRHDNMNWKTWPIFNDRDTPNGNG TNLYGAQTFFLCLEDASGLSFGVFLMNSNAMEVVLQPAPAITYRTIGGILDFYV FLGNTPEQVVQEYLELIGRPALPSYWALGFHLSRYEYGTLDNMREVVERNRA AQLPYDVQHADIDYMDERRDFTYDSVDFKGFPEFVNELHNNGQKLVIIVDPAIS NNSSSKPYGPYDRGSDMKIWNSSDGVTPLIGEVWPGQTVFPDYTNPNCAV WWTKEFELFHNQVEFDGIWIDMNEVSNFVDGGSVSGCSTNNLNNPPFTPRILD GYLFCKTLCMDAVQHWGKQYDIHNLYGYSMATAEAAKTVFPNKRSFILTRS TFAGSGKFAAHWLGDNATWDDLRSIPGVLEFNLFGIPMVGPDICGFALDTP EELCRRWMQLGAFYPFSRNHNGQGYKDQDPASFGADSLNNSRHYLNIRYT LLPYLYTLFFRAHSRGDTVARPLLHEFYEDNSTWDVHQQLWGPGLLITPVLD EGAEKVMAYVPDAVWYDYETGSQVRWRKQKVEMELPGDKIGLHLRGGYIFPT QQPNTTTLASRKNPLGLIILDENKEAKGELFDWNGETKDTVANKVYLLCEFSV TQNRLEVNISQSTYKDPNNLAFNEIKILGTEEPSNVTVKHNGVPSQTSPTVTYD SNLKVAIITDIDLLLGEAYTVEW
Research Area	Others
Source	Homo sapiens (Human)
Target Names	MGAM
Expression Region	87-954aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at




4°C for up to one week.

Tag Info C-terminal hFc1-tagged

Mol. Weight 127.2 kDa

Protein Length Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.