




Recombinant Human CD63 antigen (CD63)-VLPs

| | |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MP004950HU(A4) |
| Abbreviation | Recombinant Human CD63 protein-VLPs |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | P08962 |
| Form | Lyophilized powder |
| Storage Buffer | Lyophilized from a 0.2 μm filtered PBS, 6% Trehalose, pH 7.4. |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Sensitivity | Not Test |
| Purity | The purity information is not available for VLPs proteins. |
| Sequence | MAVEGGMKCVKFLLYVLLLAFCACAVGLIavgvgaqlvlsqtiiqgatpgsllp VVIIavgvflflvafvgccgackenyclmitfaiFLSLIMLVEVAAAIAgyvfrdk VMSEFNNFRQQMENYPKNNHTASILDRMQADFKCCGAANYTDWEKIPSMS KNRVPDSCCINVTVGCGINFNEKAIHKEGCVEKIGGWLRKNVLLVAAAALGIAF VEVLGIVFACCLVKSIRSGYEVm |
| Research Area | Cancer |
| Source | Homo sapiens (Human) |
| Target Names | CD63 |
| Expression Region | 1-238aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | C-terminal 10xHis-tagged(This tag can be tested only under denaturing conditions) |
| Mol. Weight | 27.0 kDa |
| Protein Length | Full Length |
| Image |  CSB-MP004950HU(A4) is detected by Mouse anti-6*His monoclonal antibody. |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing. |



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.