




Recombinant Human C-C chemokine receptor type 3 (CCR3)-VLPs

Product Code	CSB-MP004842HU
Abbreviation	Recombinant Human CCR3 protein-VLPs
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P51677
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered PBS, 6% Trehalose, pH 7.4.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	The purity information is not available for VLPs proteins.
Sequence	MTTSLDVTVEFTGTTSSYYDDVGLLCEKADTRALMAQFVPPPLYSLVFTVGLLGNV VVVMILIKYRRLRIMTNIYLLNLAISDLLFLVTLPFWIHYVRGHNWVFGHGMCKLL SGFYHTGLYSEIFFIILLTIDRYLAIVHAVFALRARTVTFGVITSIVTWGLAVLAALP EFIFYETEELFEETLCSALYPEDTVYSWRHFHTLRMTIFCLVLP LLVMAICYTGII KTLRLCPSKKKYKAIRLIFVIMAVFFIFWTPYNVAILLSSYQSILFGNDCERSKHL DLVMLVTEVIAYSHCCMNPVIYAFVGERFRKYLRHFFHRHLLMHLGRYIPFLPS EKLERTSSVSPSTAEPELSIVF
Research Area	Immunology
Source	Homo sapiens (Human)
Target Names	CCR3
Expression Region	1-355aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged(This tag can be tested only under denaturing conditions)
Mol. Weight	42.4 kDa
Protein Length	Full Length
Image	



CSB-MP004842HU is detected by Mouse anti-6*His monoclonal antibody.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. Aliquot for long-term storage at -80°C. Solubilize for 60 minutes at room temperature with occasional gentle mixing. Avoid vigorous shaking or vortexing.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.