





Recombinant Mouse Ly6/PLAUR domain-containing protein 3 (Lypd3) (Active)

Product Code	CSB-MP849681MO
Abbreviation	Recombinant Mouse Lypd3 protein (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91YK8
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Biological Activity	Measured by its binding ability in a functional ELISA. Immobilized Mouse Lypd3 at 2 µg/mL can bind Anti-LYPD3 recombinant antibody(CSB-RA013263MA2HU).The EC50 is 2.143-2.783 ng/mL.
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	LECYSCVQKADDGCSPHRMKTVKCGPGVDVCTEAVGAVETIHGQFSVAVRG CGSGIPGKNDRGLDLHGLLAFFQLQQCEDRCNAKLNLTLRGLNPAGNESAY EPNGAECYSCVGLSREKCGSMPPVVNCYNASGRVYKGCDFGNVTLTAANV TVSLPVRGCVQDETCTRDGVTGPGFTLSGCCQGPCNADLRNKTYFSPRIP PLVLLPPTTAAPSTRAQNSSSTTSTAAPTTTTSIKPTTAQASHTSPHEMDLEV IQEEGASLSGGAAGHGTTAGHGGAAGHQDRSNMEKYPGKGGGAQIPAKGGS
Source	Mammalian cell
Target Names	Lypd3
Expression Region	33-343aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	33.4 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity
 Measured by its binding ability in a functional ELISA. Immobilized Mouse Lypd3 at 2 µg/ml can bind Anti-LYPD3 recombinant antibody (CSB-RA013263MA2HU). The EC₅₀ is 2.143-2.783 ng/mL.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.