





Recombinant Mouse Probetacellulin (Btc), partial (Active)

| | |
|----------------------------|---|
| Product Code | CSB-MP002853MO |
| Abbreviation | Recombinant Mouse Btc protein, partial (Active) |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?. |
| Uniprot No. | Q05928 |
| Form | Lyophilized powder |
| Storage Buffer | Lyophilized from a 0.2 μm sterile filtered PBS, 6% Trehalose, pH 7.4 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Biological Activity | The ED50 as determined by the dose-dMPendent stimulation of the proliferation of NIH-3T3 cells is 0.3304-0.5261 ng/mL. |
| Purity | ≥ 95% as determined by SDS-PAGE. |
| Sequence | DGNTTRTPETNGSLCGAPGENCTGTTTPRQKVKTHFSRCPKQYKHYCIHGRCR FVVDEQTPSCICEKGYFGARCERVDLFY |
| Source | Mammalian cell |
| Target Names | Btc |
| Expression Region | 32-111aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week. |
| Tag Info | C-terminal 10xHis-tagged |
| Mol. Weight | 10.7 kDa |
| Protein Length | Partial |

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.


Activity

The ED₅₀ as determined by the dose-dependent stimulation of the proliferation of NIH-3T3 cells is 0.3304-0.5261 ng/mL.

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.