




Recombinant Human Ran-binding protein 3-like (RANBP3L)

Product Code	CSB-YP803130HU
Abbreviation	Recombinant Human RANBP3L protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q86VV4
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	MTTIPRKGSSHLPGSLHTCKLKLQEDRRQKEKS VIAQPIFVFEKGEQTFKRPAE DTLYEAAEPECNGFP TKRVRSSSFTFHITDSQSQGVRKNNVFM TSALVQSSVD IKSAEQGPVKHSHKHVIRPAILQLPQARSCAKVRKTFGHKALESCKTKEKTNNKI SEGNSYLLSENLSRARISVQLSTNQDFLGATSVGCQPNE DKCSFKSCSSNFVF GENMVERVLGTQKLTQPQLEND SYAKEKPFKSIPKFPVNFLSSRTDSIKNTSLI ESAAAFSSQPSRKCLLEKIDVITGEETEHNVLKINCKLFI FNKTTQSWIERGRGT LRLNDTASTDCGTLQSR LIMRNQGSRLILNSKLWAQM KIQRANHKNVRITATD LEDYSIKIFLIQASAQDTAYLYAAIHHRLVALQSFNKQRDVNQAESLSETAQQLN CESCDENEDDFIQVTKNGSDPSSWTHRQSVACS
Research Area	Cell Biology
Source	Yeast
Target Names	RANBP3L
Expression Region	1-465aa
Tag Info	C-terminal 6xHis-Flag-tagged
Mol. Weight	54.7 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.