



# CD19 (FMC63 Biosimilar) Recombinant Monoclonal Antibody

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>       | CSB-RA004888MB1HU                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Storage</b>            | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Uniprot No.</b>        | P15391                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Immunogen</b>          | Recombinant Human CD19 protein                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Species Reactivity</b> | Human                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Form</b>               | Liquid                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Conjugate</b>          | Non-conjugated                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Storage Buffer</b>     | 0.01M PBS,pH7.4                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Clonality</b>          | Monoclonal                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Product Type</b>       | Biosimilar Antibody                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Immunogen Species</b>  | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Lead Time</b>          | In stock                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Target Names</b>       | CD19                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Notes</b>              | <p>Validation Status<br/>Application-specific performance (e.g., in flow cytometry, ELISA, IHC or other assay formats) has not yet been experimentally verified by CUSABIO. Users are advised to determine the optimal working conditions empirically in their own assay systems.</p> <p>Guaranteed Quality</p> <p>① Antibody purity &gt; 95% tested by SDS-PAGE.</p> <p>② Endotoxin level &lt; 0.1EU/ug tested by LAL method.</p> |