



IL12&IL23 (Ustekinumab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA011587MB1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29459&P29460
Immunogen	Recombinant Human IL12&IL23 protein
Species Reactivity	Human
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS,pH7.4
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	In stock
Target Names	IL12&IL23

Notes

Validation Status

Application-specific performance (e.g., in flow cytometry, ELISA, IHC or other assay formats) has not yet been experimentally verified by CUSABIO. Users are advised to determine the optimal working conditions empirically in their own assay systems.

Guaranteed Quality

- ① Antibody purity > 95% tested by SDS-PAGE.
- ② Endotoxin level < 0.1EU/ug tested by LAL method.