



SLC39A6 (Ladiratuzumab vedotin Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA621669MB1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13433
Immunogen	Recombinant Human SLC39A6 protein
Species Reactivity	Human
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS,pH7.4
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	In stock
Target Names	SLC39A6

Notes

Validation Status

Application-specific performance (e.g., in flow cytometry, ELISA, IHC or other assay formats) has not yet been experimentally verified by CUSABIO. Users are advised to determine the optimal working conditions empirically in their own assay systems.

Guaranteed Quality

① Antibody purity > 95% tested by SDS-PAGE.

② Endotoxin level < 0.1EU/ug tested by LAL method.