



TNFSF11 (Denosumab-dssb Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA023986MB5HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14788
Immunogen	Recombinant Human TNFSF11 protein
Species Reactivity	Human
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS,pH7.4
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	3-4 weeks
Target Names	TNFSF11

Notes

Validation Status

Application-specific performance (e.g., in flow cytometry, ELISA, IHC or other assay formats) has not yet been experimentally verified by CUSABIO. Users are advised to determine the optimal working conditions empirically in their own assay systems.

Guaranteed Quality

- ① Antibody purity > 95% tested by SDS-PAGE.
- ② Endotoxin level < 0.1EU/ug tested by LAL method.