



# TNFRSF4 (INCAGN-1949 Biosimilar) Recombinant Monoclonal Antibody

|                           |   |
|---------------------------|---|
| <b>Product Code</b>       | CSB-RA023981MB10HU  |
| <b>Storage</b>            | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| <b>Uniprot No.</b>        | P43489  |
| <b>Immunogen</b>          | Recombinant Human TNFRSF4 protein                             |
| <b>Species Reactivity</b> | Human   |
| <b>Form</b>               | Liquid  |
| <b>Conjugate</b>          | Non-conjugated  |
| <b>Storage Buffer</b>     | 0.01M PBS,pH7.4   |
| <b>Clonality</b>          | Monoclonal  |
| <b>Product Type</b>       | Biosimilar Antibody   |
| <b>Immunogen Species</b>  | Homo sapiens (Human)  |
| <b>Lead Time</b>          | 3-4 weeks   |
| <b>Target Names</b>       | TNFRSF4   |

## Notes

### Validation Status

Application-specific performance (e.g., in flow cytometry, ELISA, IHC or other assay formats) has not yet been experimentally verified by CUSABIO. Users are advised to determine the optimal working conditions empirically in their own assay systems.

### Guaranteed Quality

- ① Antibody purity > 95% tested by SDS-PAGE.
- ② Endotoxin level < 0.1EU/ug tested by LAL method.