



CD40 (Ravagalimab Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA004936MB18HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25942
Immunogen	Recombinant Human CD40 protein
Species Reactivity	Human
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS,pH7.4
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	3-4 weeks
Target Names	CD40
Notes	<p>Validation Status Application-specific performance (e.g., in flow cytometry, ELISA, IHC or other assay formats) has not yet been experimentally verified by CUSABIO. Users are advised to determine the optimal working conditions empirically in their own assay systems.</p> <p>Guaranteed Quality</p> <p>① Antibody purity > 95% tested by SDS-PAGE.</p> <p>② Endotoxin level < 0.1EU/ug tested by LAL method.</p>