



MS4A1 (Rituximab-arrx Biosimilar) Recombinant Monoclonal Antibody

Product Code	CSB-RA015007MB19HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11836
Immunogen	Recombinant Human MS4A1 protein
Species Reactivity	Human
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	0.01M PBS,pH7.4
Clonality	Monoclonal
Product Type	Biosimilar Antibody
Immunogen Species	Homo sapiens (Human)
Lead Time	3-4 weeks
Target Names	MS4A1
Notes	<p>Validation Status Application-specific performance (e.g., in flow cytometry, ELISA, IHC or other assay formats) has not yet been experimentally verified by CUSABIO. Users are advised to determine the optimal working conditions empirically in their own assay systems.</p> <p>Guaranteed Quality</p> <p>① Antibody purity > 95% tested by SDS-PAGE.</p> <p>② Endotoxin level < 0.1EU/ug tested by LAL method.</p>