




# Recombinant Human C-C chemokine receptor type 7 (CCR7)-VLPs

<b>Product Code</b>	CSB-MP004846HU
<b>Abbreviation</b>	Recombinant Human CCR7 protein-VLPs
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P32248
<b>Form</b>	Lyophilized powder
<b>Storage Buffer</b>	Lyophilized from a 0.2 μm filtered PBS, 6% Trehalose, pH 7.4.
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	VLPs
<b>Sequence</b>	QDEVTDYIGDNTTVDYTLFESLCSKKDVRNFKAWFLPIMYSIICFVGLLGNGL VVLTYIYFKRLKTMTDTYLLNLAVADILFLLTLPFWAYSAAKSWVFGVHFCKLIF AIYKMSFFSGMLLLLCISIDRYVAIVQAVSAHRHRARVLLISKLSVGIWILATVLS IPELLYSDLQRSSESEQAMRCSLITEHVEAFITIQVQMVGFLVPLLAMSFCYLVII RTLLQARNFERNKAIKVIIAVVVVVFVFLPYNGVVLAQTVANFNITSSTCELSKQ LNIAYDVTYSLACVRCVNPFLYAFIGVKFRNDLFLKFKDLGCLSQEQLRQWSS CRHIRRSSMSVEAETTTTFSP
<b>Research Area</b>	Others
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CCR7
<b>Expression Region</b>	25-378aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	C-terminal 10xHis-tagged
<b>Mol. Weight</b>	41.8 kDa
<b>Protein Length</b>	Full Length of Mature Protein

Image



CSB-MP004846HU is detected by Mouse anti-6\*His monoclonal antibody.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.