




Recombinant Human Growth/differentiation factor 8 (MSTN)

Product Code	CSB-MP015057HU(A4)d7
Abbreviation	Recombinant Human MSTN protein
Uniprot No.	O14793
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.61.
Immunogen Species	Homo sapiens (Human)
Purity	≥ 95% as determined by SDS-PAGE.
Sequence	NENSEQKENVEKEGLCNACTWRQNTKSSRIEAIKIQILSKLRLETAPNISKD VIR QLLPKAPPLRELIDQYDVQRDDSSDGSLEDDDYHATTETIITMPTESDFLMQVD GKPKCCFFKFSSKIQYNKVVKAQLWIYLRPVETPTTVFVQILRLIKPMKDGTRYT GIRSLKLDMNPGTGIWQSIDVKTVLQNW LKQPESNLGIEIKALDENGHD LAVTF PGPGEDGLNPFLEVKVTDTPKRSRRDFGLDCDEHSTESRCCRYPLTVDFEAF GWDWIIAPKRYKANYCSGECEVFVLQKYPHTHLVHQANPRGSAGPCCTPTKM SPINMLYFNGKEQIIYGKIPAMVVDRCGCS
Research Area	Signal Transduction
Source	Mammalian cell
Target Names	MSTN
Expression Region	24-375aa
Tag Info	C-terminal 10xHis-tagged
Mol. Weight	42.1 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.