




# Recombinant Human Galectin-9 (LGALS9)(W309A)

<b>Product Code</b>	CSB-MP012895HU(F1)(M)d7
<b>Abbreviation</b>	Recombinant Human LGALS9 protein (W309A)
<b>Uniprot No.</b>	O00182
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.67.
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 95% as determined by SDS-PAGE.
<b>Sequence</b>	MAFSGSQAPYLSPAVPFSGTIQGGLQDGLQITVNGTVLSSSGTRFAVNFQTGF SGNDIAFHFNPRFEDGGYVVCNTRQNGSWGPEERKTHMPFQKGMPPDLCFL VQSSDFKVMVNGILFVQYFHRVPFHRVDTISVNGSVQLSYISFQNPRTVPVQP AFSTVPFSQPVCFPPRPRGRRRQKPPGVWPANPAPITQTVIHTVQSAPGQMFS TPAIPPMYPHPAYPMPFITLGGLYPSKSILLSGTVLPSAQRFHINLCSGNHIA FHLNPRFDENAVVRNTQIDNSWGSEERSLPRKMPFVRGQSFSVAILCEAHCLK VAVDGQHLFEYYHRLRNLPTINRLEVGGDIQLTHVQT
<b>Research Area</b>	Neuroscience
<b>Source</b>	Mammalian cell
<b>Target Names</b>	LGALS9
<b>Expression Region</b>	1-355aa(W309A)
<b>Tag Info</b>	C-terminal 10xHis-tagged
<b>Mol. Weight</b>	42.1 kDa
<b>Protein Length</b>	Full Length

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.